
Information on handling, possibilities of use and any restrictions for the product group:

Trapeziometacarpal joint replacement

ELIS – Cementless stem, T	ELIS – Cementless stem, T/II (Ti+HA)
Size range: 1 – 5	Size range: 1 – 5
ELIS – Cementless stem (CC)	ELIS – Straight neck, T/II
Size range: 1 – 5+	Size range: L12 – 18
ELIS – Valgus neck 15°, T/II	ELIS – Valgus neck 15°/Ante-L/R
Size range: L12 – 18	Size range: R/L, L12 – 18
ELIS – DM – Straight neck	ELIS – DM – Valgus neck 15°
Size range: L12 – 18	Size range: L12 – 18
ELIS – DM – Valgus neck 15°/Ante-L/R	ELIS – Cemented cup with offset, T/II
Size range: R/L, L12 – 18	ELIS – Cementless cup, T/III
ELIS – Sphere – Cemented cup	ELIS – Cementless cup, T/III
ELIS – Sphere – Cementless cup (CC)	ELIS – Insert, T/III
Size range: D9 – 10	ELIS – Sphere, DM – Cementless cup (CC)
ELIS – Sphere, DM – Cementless cup (CC)	ELIS – Konos, DM – Cementless cup (CC)
Size range: D9 – 10	Size range: D9 – 10
ELIS – Trapezium replacement, TR	ELIS – Trapezium insert, TR
Size range: 14/15 – 16/17	Size range: 9/5 – 11/5

Intended purpose

Implants are intended as a replacement for the trapeziometacarpal joint.

Product description

The trapeziometacarpal joint replacement in the stem, neck and cup assembly is intended for primoimplantation in cases of damage to the articulating parts of the trapezium and metacarpus that can no longer be solved by conservative treatment. In special cases, the implant can also be used as a revision implant in cases of failure of the primary implant (revision surgery is decided by the physician on the basis of a preoperative examination). In this set, the trapeziometacarpal joint replacement always consists of three parts – the metacarpal stem, the neck with the head and the cementless cup placed in the trapezium (cup + PE insert) or the cemented cup.

The trapeziometacarpal joint replacement in the assembly of stem, neck and trapezium replacement is intended for revision surgeries when it is no longer possible to implant the cup safely due to defects in the trapezium. In addition, it is intended to solve severe Eaton-Littler IV arthrosis when trapezium removal is necessary.

In this set-up, the trapeziometacarpal joint replacement always consists of three parts – the metacarpal stem, the neck with a 5 (mm) diameter head and the trapezium replacement (body + PE insert). The replacement cannot be used as a partial replacement, only as a full replacement.

ELIS – Cementless stem, T, ELIS – Cementless stem, T/II (Ti+HA) and ELIS – Cementless stem (CC)
The cementless stem is designed to follow the shape and dimensions of the metacarpal cavity. Body of the stem is slightly curved and its cross-section is triangular with rounded edges. For the insertion and fixation of the neck, it has a conical hole with an oval recess at the top to ensure the anti-rotational stability of the neck. An M3 thread is below the bottom of the cone for the extractor. The stem is made of a wrought titanium 6-aluminum 4-vanadium alloy (ISO 5832-3). The outer surface in the upper part is coated with plasma sprayed unalloyed titanium (ISO 5832-2), double layer unalloyed + hydroxyapatite Ti + HA (ISO 5832-2 + ISO 13779-2) or composite coating unalloyed titanium with hydroxyapatite Ti with HA (ISO 5832-2 with ISO 13779-2). The lower end of the stem in about 1/3 of the length is finely blasted. Stems are designed for assembly with all types of necks.

ELIS – Straight neck, T/II, ELIS – Valgus neck 15°, T/II and ELIS – Valgus neck 15°/Ante-L/R

The neck (ELIS) is available in straight and 15° valgus versions. In addition, the neck (ELIS) in the 15° valgus version is extended with an anteverision, right and left (R/L). They differ from each other by the position of the axis of the neck against the axis of the connecting conical hole in the stem. Each version has four lengths (12, 14, 16 and 18 mm). All necks have 5 mm diameter head and are made of wrought cobalt-chromium-molybdenum alloy (ISO 5832-12). The surface is brushed except for the articulation surface of the head, which is polished.

ELIS – DM – Straight neck, ELIS – DM – Valgus neck 15° and ELIS – DM – Valgus neck 15°/Ante-L/R

Dual mobility necks are in straight, 15° valgus and 15° valgus in anteverision, right and left (R/L) versions. They differ in the position of the neck axis as the above mentioned ELIS necks. Each type of neck has four lengths (12, 14, 16 and 18 mm). The 5 mm diameter articulating head is fitted with a PE insert which is an integral part of the neck. The outer spherical surface of the PE insert serves to articulate with the inner spherical surface of the ELIS – DM, Sphere or Konos cup. The neck is made of wrought cobalt-chromium-molybdenum alloy (ISO 5832-12). The surface of the neck is brushed except for the articulation surface of the head, which is polished. The PE neck liner is made of crosslinked Ultra-high-molecular-weight polyethylene (ISO 5834-2).

ELIS – Cementless cup, ELIS – Sphere – Cementless cup (CC) and ELIS – Insert, T/III
Both cups consist of two separate implants (cup and PE insert). They are designed for fixation in the first metacarpal joint (in trapezium) without the use of bone cement.

The **cementless cup** has spherical shape with a diameter of 9 or 10 mm. The outer surface has anti-rotation elements on the upper part of outer layer and a porous unalloyed titanium and hydroxyapatite Ti + HA double coating (ISO 5832-2 + ISO 13779-2) (T/III) or a porous unalloyed titanium Ti composite coating with hydroxyapatite HA (CC) (ISO 5832-2, ISO 13779-2) is applied to the entire surface. The inner spherical surface is designed to accommodate the PE insert and has elements to secure the insert against falling out. The cup is designed for insertion of guide wire and has a 2 mm diameter central hole at the bottom. The cup is made of a wrought titanium 6-aluminum 4 vanadium alloy (ISO 5832-3). The **PE insert** has a spherical outer surface with a diameter corresponding to the cup dimension. Locking elements designed to connect with the cup are at the top of insert. The inside of PE insert is designed to articulate with a 5 mm diameter head of neck. The PE insert is made of cross-linked ultra-high molecular weight polyethylene (ISO 5834-2).

ELIS – Sphere, DM – Cementless cup (CC) and ELIS – Konos, DM – Cementless cup (CC)

The cementless DM cup (Sphere, Konos) consists of a cup (Sphere, Konos) into which the DM neck with a PE insert is inserted during the surgery. It is designed to be fixed in the trapezium without the need of cement. The cementless ELIS – Sphere cup is spherical and Konos is conical. Both cups come in 9 and 10 mm sizes. The internal articulating spherical surface is designed for the PE neck insert, which is an integral part of the neck (ELIS – DM) with a 7.7 mm diameter head. The cup (ELIS-DM) does not have a hole for insertion of guide wire. The cup is made of wrought cobalt-chromium-molybdenum alloy (ISO 5832-12). The outer surface is coated

with an osteointegrating porous composite coating of unalloyed titanium and hydroxyapatite Ti + HA (ISO 5832-2, ISO 13779-2).

ELIS – Cemented cup with offset, T/II and ELIS – Sphere – Cemented cup

Cups are designed as cemented. They are designed for the II assembly. The outer surface of the cup is a hemispherical surface with a diameter of 10 mm. The inner articulation surface is spherical with 5 mm diameter. There are stabilising grooves on the perimeter of the cup designed for better bonding with the cement. They are made of cross-linked ultra-high molecular weight polyethylene (ISO 5834-2). In addition, the cup (ELIS) has a diagnostic element on its outer surface for easier identification of the position on the X-ray image. The diagnostic wire is made of wrought stainless steel (ISO 5832-1).

ELIS – Trapezium repl., TR and ELIS – Trapezium insert, TR

The trapezium replacement consists of two separate parts (trapezium replacement body and PE insert) and is designed for complete trapezium replacement in a total trapeziometacarpal joint replacement assembly.

The cementless **body of the trapezium replacement** has a barrel-shaped elliptical cross-section. It is manufactured in three sizes 14/15, 15/16 and 16/17 mm and all of them have 10 mm height. The outer surface has a groove on one side for stabilisation by the radial part of the flexor carpi radialis and two transverse holes for further ancillary fixation. The inner part of the body has tapered hole for the PE insert.

The **PE insert** has a tapered outer surface with dimensions corresponding to the dimensions of the body of the trapezium replacement. The inner spherical surface is designed for a neck with 5 mm diameter head. The body of the replacement is made of wrought cobalt-chromium-molybdenum alloy (ISO 5832-12). The outer surface is polished. The PE insert is made of cross-linked ultra-high molecular weight polyethylene (ISO 5834-2).

Intended Patient Population (Target Group)

- Age - 18+ adult population, male, female
- Weight - without limitation.
- Health - the patient must be mentally capable of understanding the significance of the surgery and adhering to the post-operative regimen. The patient must be aware of the limitations caused by the implant.
- Nationality - not relevant.

Indications

- Rhizarthrosis primary and secondary, resistant to conservative treatment.
- Relative indication is rheumatoid arthritis (not valid for trapezium replacement).
- X-ray criterion Eaton-Littler II-IV.
- Unfixed 2 deformity of the thumb at the metacarpal joint.

Contraindications

- Young, very active patient (incomplete growth).
- X-ray Eaton-Littler IV (not valid for trapezium replacement).
- A relative contraindication is a rigid Z deformity of the thumb at the metacarpal joint (to be treated bilaterally or simultaneously with implant of the trapeziometacarpal joint).
- Small trapezium (not valid for trapezium replacement).
- Trapezium dysplasia (Kapandji) >140° (not valid for trapezium replacement).
- Severe osteoporosis.
- Allergies (Ti, Al, V, Co, Cr, Mo, UHMWPE, hydroxyapatite, bone cement).
- Inflammatory diseases specific and non-specific.
- Tumours in the relevant localisation.
- Patient's inability to cooperate after surgery.

Clinical benefits that can be expected

- Dual Mobility – statistically significant improvement in outcomes and significant improvement in joint stability (due to the own design of the medical device)
- Reduction of pain – VAS score compared to preoperative examination – at least 45% lower value.
- QuickDASH score reduced by at least 45%.
- Increase in range of motion (opposition according to Kapandji) – at least 1 point more on the scale.
- Key Pinch – at least 35% more force.
- Increased Hand Grip Strength – at least 10% more force.
- Osteointegration about 30 % better with composite coating than with Ti+HAp coating.

Reference to the Summary of Safety and Clinical Performance (SSCP): available in the European Medical Device Database (Eudamed) where it is linked to the Basic UDI-DI – 85926020700FZ. URL of the Eudamed public website: <https://ec.europa.eu/tools/eudamed> or on the website www.beznoska.cz.

Properties of the product in terms of effectiveness:

Implants fully replace the trapeziometacarpal joint.

Residual risks:

Danger of failure of the medical device if the user does not follow the requirements for use set out in the instruction for use.

Qualifications:

- Doctor – Orthopedist, Traumatologist, Surgeon.
- Experience with trapeziometacarpal joint implantation.


Bone cement:

Different types of bone cement are used in the cemented version of trapeziometacarpal joint replacement. BEZNOSKA company supplies two types of bone cement: BEZCEM – Surgical Cement, BEZGEM – Surgical Cement with Gentamicin. It is possible that in hospitals, in addition to the bone cements supplied by BEZNOSKA, other bone cements may be encountered. The condition for use is that the bone cement meets the basic standard ISO 5833.

Materials used:

(Standard / Number of standard / Name of standard)
ISO 5832-1 – Implants for surgery – Metallic materials – Part 1: Wrought stainless steel
ISO 5832-2 – Implants for surgery – Metallic materials – Part 2: Unalloyed titanium
ISO 5832-3 – Implants for surgery – Metallic materials – Part 3: Wrought titanium 6-aluminum 4-vanadium alloy
ISO 5832-12 – Implants for surgery – Metallic materials – Part 12: Wrought cobalt-chromium-molybdenum alloy
ISO 5834-2 – Implants for surgery – Ultra-high-molecular-weight polyethylene – Part 2: Moulded forms
ISO 13779-2 – Implants for surgery – Hydroxyapatite – Part 2: Thermally sprayed coatings of hydroxyapatite

Combinations

Warnings and cautions

Before surgery

Patient

- The patient consents to the surgery and associated risks – see Adverse Complications.
- The patient must be mentally able to understand the significance of the surgery and comply with the postoperative regime.
- The patient must be aware of the limitations caused by the implant.

Departments performing the operation

- Medical qualifications – orthopedist, surgeon, traumatologist.
- Ensure highly sterile surgical conditions.
- Check the sterility of the packaging (packaging damage and expiration date).
- Instruments for implantation must be complete and functional.
- Use the instrumentation supplied by the manufacturer, do not use others.
- The implant is designed for single use only.
- Ensure that damaged, unsterile or reoperated implant is not used.
- In the event of damage to the outer packaging, contact the manufacturer at the address listed at the end of this instruction for use.
- In case of damage to the inner packaging, dispose the implant as hazardous waste – see Disposal of the product.
- Do not use the implant with patients with documented metal allergy or intolerance (Ti, Al, V, Co, Cr, Mo, UHMWPE, HAp, bone cement).
- Patients who receive a trapeziometacarpal joint replacement must be advised that the lifespan of the implant is dependent on their weight and their activity level.
- All necessary implant components must be available.

During the surgery

- During the surgery it is necessary to use the instrumentation supplied by the manufacturer and to follow the surgical procedure recommended by the manufacturer.
- It is necessary to use surgical gloves when handling the implant.
- Incorrect resection or incorrect positioning of components must be avoided.
- When handling the product, the polished articulation surfaces, surfaces coated with Ti or Ti + HA plasma coating and the threaded surfaces of the stem must be protected. Damaged articulation surfaces adversely affect the life of the implant.
- Damaged instruments and devices must not be used during the surgery.
- The neck implants must not be used for trial bending.
- When using a cemented cup assembly, remove all loose bone cement; pieces of bone cement may penetrate between the friction surfaces and cause abnormal wear of the implant.
- Always insert the implant into a precisely machined (verified by trial components), clean and dry bed of the bone.
- Articulation surfaces must be perfectly clean during the trial and final articulation.
- Clean and dry the articulation surfaces thoroughly before inserting the final neck.

After the surgery

Patient

- The patient must be aware of the limitations caused by the implant (In the first three months, follow a strict postoperative regimen – limit exercise and avoid any vibrations).
 - Patient must be warned that his new implant can only be subjected to a limited load until the bones are fully healed.
 - In case of unexpected changes (any change in health, i.e. pregnancy) that may be related to the implant, it is recommended to see a specialist.
 - Based on the available information, MRI examinations with a static magnetic field up to 3T can be performed in patients with implants from BEZNOSKA, s.r.o., but not earlier than 6 weeks after implantation and without signs of loosening (see Declaration on the possibility of MRI examinations + methodological sheet for MRI examinations of patients with metal implants + contraindications and risks of MRI examinations – Information can be found on the website <https://www.beznoska.cz/pro-pacienty/nejcastejsi-otazky>)
- Department performing the surgery
- The load-bearing capacity of the implant cannot be compared to the healthy bone.
 - Timely preventive checks of the implant.

Surgery procedure

The trapeziometacarpal joint has its own detailed surgical procedure manual, where the complete implantation is described.

Adverse Complications

During surgery

- Damage to vascular and nervous structures.
- Iatrogenic bone damage, sometimes leading to fracture.
- Incorrect position of components, incorrect resection.

After surgery

- Cardiovascular disorders such as VTE, deep phlebothrombosis, MI, postoperative hematoma.
- Wound healing disorders, infection.
- Loosening of the entire total replacement or its component, deformation or breakage of the component.

- Adverse reactions of the organism to implant.

Disposal of the product

After use, the implant is listed as „hazardous waste“ for the healthcare. Disposal according to the applicable legislation of the country where the implant is being disposed.

Manufacturer's regulation

The implant cannot be sterilised and used again, because of following risks:

- Damage to the polished articulating, connecting and anchoring surfaces.
- Damage or contamination of the porous coating.

Risks of reuse from the patient's perspective

- Implant failure due to the damage to articulating, connecting and anchoring surfaces, or improper sterilization – risk of early reoperation.
- Risk of infection in the patient – risk of difficult to treat infection and subsequent loosening.
- Change of implant size or type, or incompleteness – risk of re-operation or early reoperation.

The implant information is provided in the implant card in accordance with Article 18 of MDR 2017/745. The card is attached to the medical device (further information on the implant card at <https://www.beznoska.cz/pro-pacienty/nejcastejsi-otazky/> in the „karta implantátu“ section).

Reporting of suspected adverse reactions/adverse events:

Patients and healthcare professionals should be instructed on how to proceed when reporting adverse events and suspected adverse events or any problems related to the medical device. Any serious adverse event occurring in relation to the medical device should be reported to the manufacturer (at the address given at the end of this instructions) and to the competent authority of the Member State where the user or patient is established.

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
	Manufacturer		Date of manufacture
	Use by date		Temperature limit
	Humidity limitation		Protect from heat and radioactive sources
	Do not re-use		Do not resterilize
	Medical device		Unique device identifier
	Catalogue number		Batch code
	Quantity		Sterilized using ethylene oxide
	Cemented		Cementless
	Double sterile barrier system		Material
	Caution		Date
	Patient identification		Health care centre or doctor
	Contains hazardous substances		
	Operator's manual; operating instructions		
	Refer to instruction manual/booklet		
	Do not use if package is damaged		
	Patient information website		

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ELIS – Cementless cup, ELIS – Sphere – Cementless cup (CC) og ELIS – Insert, T/II

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ELIS – Sphere, DM – Cementless cup (CC) og ELIS – Konos, DM – Cementless cup (CC)

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ELIS – Cemented cup with offset, T/II og ELIS – Sphere – Cemented cup

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
Anvendte materialer:

(Standard / Standardens nummer / Standardens titel)
 ISO 5832-1 – Implants for surgery – Metallic materials – Part 1: Wrought stainless steel
 ISO 5832-2 – Implants for surgery – Metallic materials – Part 2: Unalloyed titanium
 ISO 5832-3 – Implants for surgery – Metallic materials – Part 3: Wrought titanium 6-aluminium 4-vanadium alloy
 ISO 5832-12 – Implants for surgery – Metallic materials – Part 12: Wrought cobalt-chromium-molybdenum alloy
 ISO 5834-2 – Implants for surgery – Ultra-high-molecular-weight polyethylene – Part 2: Moulded forms
 ISO 13779-2 – Implants for surgery – Hydroxyapatite – Part 2: Thermally sprayed coatings of hydroxyapatite

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