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

Information on handling, possibilities of use and any restrictions for the product group:

TOTAL TRAPEZIOMETACARPAL JOINT REPLACEMENT

Identification

ELIS – Cementless stem, T	ELIS – Cementless stem, T/II (Ti+HA)
Size range: 1 – 5	Size range: 1 – 5
ELIS – Cementless stem (CC)	ELIS – Straight neck, T/II
Size range: 1 – 5	Size range: L12 – 18
ELIS – Valgus neck, T/II	ELIS – Valgus neck 15°/ANTE-R, T/II
Size range: L12 – 18	Size range: L12 – 18, L/R
ELIS – DM – Straight neck, T/III	ELIS – Cementless cup, T/III
Size range: L12 – 18, L/R, straight/valgus	Size range: D9 – D10
ELIS – Sphere – Cementless cup (CC)	ELIS – Sphere, DM – Cementless cup (CC)
Size range: D9 – D10	Size range: D9 – D10
ELIS – Konos, DM – Cementless cup (CC)	ELIS – Insert, T/III, for head D5
Size range: D9 – D10	
ELIS – Cup with offset, T/II, 5/10	ELIS – Sphere – Cemented cup, 5/10
ELIS – Trapezium replacement, TR	ELIS – Trapezium insert, TR
Size range: 14/15 – 16/17	Size range: 9/5 – 11/5
ELIS – Cemented cup with offset, T/II	

Information about sterility

- The product is supplied by the manufacturer in a sterile pack
- Sterility of the implant in the protective sterile packaging is guaranteed for 5 years. Expiry dates are marked clearly on the protective sterile packaging
- Not to be used if the sterile packaging is damaged. The product cannot be re-sterilised

Intended use

ELIS – Total Trapeziometacarpal Joint Replacement is designed for primary replacement of the trapeziometacarpal joint in the stem, neck, cup assembly.

Product description

ELIS – Cementless stem – type T, T/II (Ti+HA) and ELIS – Cementless stem (CC)

The cementless TEP stem is designed to match the shape and dimensions of the metacarpal cavity. The body of the stem is slightly curved and its cross-section is triangular with rounded edges. For the insertion and fixation of the neck, there is a conical hole at the top with an oval recess to ensure the anti-rotational stability of the neck. Below the bottom of the cone is an M3 thread for the extractor. The stem is made of titanium alloy Ti6Al4V (ISO 5832-3). The outer surface is plasma coated on the top with Ti (ELIS type T), Ti + HA (ELIS type T/II), (ISO 5832-2 + ISO 13779-2) or a composite layer of Ti with HA (ELIS type CC), (ISO 5832-2 with ISO 13779-2). The lower end of the stem in approx. 1/3 of the length is finely blasted. The stems are designed for assembly with all types of necks.

ELIS – Neck – type T/II and ELIS – Neck straight, valg., valg./ante

The neck (ELIS) is available in straight and 15° valgus versions. In addition, the neck (ELIS) in the 15° valgus version is extended with an anteversion, right and left versions. They differ from each other by the position of the axis of the neck against the axis of the connecting conical hole in the stem. Each version has four lengths (12, 14, 16 and 18 mm). The diameter of the head is 5 mm. All types and designs are made of forged cobalt-chrome-molybdenum alloy (ISO 5832-12). The surface is brushed except for the articulation surface of the head, which is polished.

ELIS – DM – Neck

The neck is straight, 15° valgus and 15° valgus in anteversion, right and left. They differ from each other in the position of the neck axis as the above mentioned ELIS necks. Each version also has four lengths (12, 14, 16 and 18 mm). The 5 mm diameter articulating head is fitted with a PE insert, which is an integral part of the neck. The outer spherical surface of the PE insert serves to articulate with the inner spherical surface of the ELIS – DM – Sphere or Konos cup. The neck is made of forged cobalt-chrome-molybdenum alloy (ISO 5832-12). The surface of the neck is brushed except for the articulation surface of the head, which is polished. The PE neck liner is made of cross-linked ultra-high molecular weight polyethylene (ISO 5834-2).

ELIS – Cementless cup – type T/II and ELIS – Sphere – Cup

Both cups consist of two separate MDs (cup casing and PE liner). They are designed for fixation in the first metacarpal bone (os trapezium) without the use of bone cement.

The **cementless cup** is spherical with a diameter of 9 or 10 mm. The outer surface has anti-rotation elements around the perimeter at the top (at the largest diameter) and a porous Ti+HA double layer (ISO 5832-2 + ISO 13779-2), (MCP) or a porous composite layer of Ti with HA (ELIS all types), (Ti: ISO 5832-2, HA: ISO 13779-2) is applied to the entire surface. The inner spherical surface is designed to accommodate the PE insert and is provided with elements to secure the insert against falling out. The cup is adapted for insertion along the K-wire and therefore has a central hole of 2 mm diameter in the bottom. The cup is made of titanium alloy Ti6Al4V (ISO 5832-3).

The **PE insert** has a spherical outer surface with a diameter corresponding to the cup dimension. At the top there are locking elements designed to connect with the cup. The inner spherical surface of the insert is designed to articulate with a neck head of 5 (mm) diameter. The PE liner is made of cross-linked ultra-high molecular weight polyethylene (ISO 5834-2).

ELIS – Sphere or Konos, DM – Cup

The cementless (Sphere, Konos) DM consists of a cup (Sphere, Konos) into which the DM neck with a PE insert is inserted during the actual hammering. It is designed to be fixed in the trapezoid without cement.

Cup ELIS – Sphere, Konos – DM

ELIS cementless cup is spherical and Konos is conical. Both cups have a size of 9 or 10 mm. The inner articulating spherical surface is designed to accommodate the PE neck insert, which is an integral part of the neck (ELIS – DM) with a 7 mm diameter head. The cup (ELIS – DM) does not have a hole for insertion along the K-wire. The cup is made of forged cobalt-chromium-molybdenum alloy (ISO 5832-12). The outer surface is coated with an osseointegrating porous composite layer of Ti with HA (Ti: ISO 5832-2, HA: ISO 13779-2).

ELIS – Cemented cup – type T/II and ELIS – Sphere – Cemented cup

The cups are designed as cemented. They are designed for assembly (II). The outer surface of the cup is a 10 mm diameter semi-ball surface. The inner articulation surface is spherical with a diameter of 5 mm. There are stabilising grooves on the perimeter of the cup, intended for bonding with the cement. They are made of cross-linked ultra-high molecular weight polyethylene (ISO 5834-2). In addition, the cup (ELIS) has a diagnostic element on its outer surface for easier identification of the position on the X-ray image. The diagnostic wire is made of stainless steel (ISO 5832-1).

ELIS – Trapezium replacement – type TR

The trapezium replacement assembly is made of two separate parts (trapezium body and PE insert) and is designed for complete trapezium replacement in a total trapeziometacarpal joint replacement assembly.

The **cementless trapezium body** has a barrel-shaped elliptical cross-section. It is manufactured in three sizes of 14/15, 15/16 and 16/17 mm and a height of 10 mm. The outer surface is provided with a groove on one side, intended for stabilisation by the radial part of the flexor carpi radialis, and two transverse holes for further ancillary fixation. The inner part of the body is provided with a tapered hole for the PE insert.

The **PE insert** has a tapered outer surface with dimensions corresponding to the dimensions of the hole in the body of the prosthesis. The inner spherical surface is designed to articulate with a neck head having a diameter of 5 (mm). The body of the restoration is made of forged cobalt-chromium-molybdenum alloy (ISO 5832-12). The outer surface is polished. The PE liner is made of cross-linked ultra-high molecular weight polyethylene (ISO 5834-2).

Indications

- Adult population 18+, female, male
- Cases of primary and secondary rizarthritis, resistant to conservative treatment, relative indication is rheumatoid arthritis (not applicable for trapezium replacement), Eaton-Littler II-IV radiographic criteria, Unfixed Z deformity of the thumb at the MCP joint, Trapezium height (not applicable for trapezium replacement) > 7 mm


Contraindications

Young, very active patient, X-ray Eaton-Littler IV (not valid for trapezium replacement), relative contraindication is rigid Z deformity of the thumb in MCP joint (to be solved bilaterally or simultaneously with TEP of MCP joint, small trapezium (not valid for trapezium replacement), trapezoid dysplasia (Kapandji >140°) (not applicable for trapezium replacement), severe osteoporosis, documented allergy or intolerance to metals (titanium alloy ISO 5832-3, titanium ISO 5832-2), UHMWPE (ISO 5834-2) and bone cement, inflammatory diseases specific and non-specific

Combinations

TMCJ	Stem type T and T/II	1	2	3	4	5			
ELIS	Stem type (CC)	1	1+	2	2+	3	4	5	5+

TMCJ		ELIS				
Neck – type T/II		Valgus neck 15°, ante		Neck DM		
				Straight	Valgus	Ante
-Straight-	-Valgus-	-R (Right) – / -L (Left) –	-	-Straight-	-Valgus-	-R (Right) – / -L (Left) –
5/12 (S)	5/12 (S)	5/12		5/12	5/12	5/12
5/14 (M)	5/14 (M)	5/14		5/14	5/14	5/14
5/16 (L)	5/16 (L)	5/16		5/16	5/16	5/16
5/18 (XL)	5/18 (XL)	5/18		5/18	5/18	5/18

Caution and warnings

Before surgery

Patient

- The patient consents to the surgery and associated risks
- The patient must be mentally capable of understanding the significance of the operation and of complying with the post-operative regimen
- The patient must be aware of the limitations caused by the implant

Workplace performing surgery

- Ensures highly sterile operating conditions
- Checks packaging for sterility (packaging damage and expiration time)
- Instruments to be implanted must be complete and functional
- Implant is intended for single use only
- Ensures that no damaged, unsterile or reoperated implant is used
- Implant should not be used in patients with documented allergies or intolerances to metals (titanium alloy/ISO 5832-3, titanium ISO 5832-2), UHMWPE (ISO 5834-2) and bone cement
- Patients who receive a trapeziometacarpal joint replacement must be advised that the longevity of the implant is dependent on their weight and activity level
- All necessary implant components must be available

During surgery

- During application, it is necessary to use the instrumentation supplied by the manufacturer and to follow the operating procedure recommended by the manufacturer

It is necessary to use surgical gloves during manipulation

- Incorrect resection or incorrect positioning of components must be avoided
- When handling the product, the polished articulating surfaces, the Ti or Ti+HA plasma-coated surfaces and the threaded surfaces of the stem must be protected. Damaged articulation surfaces adversely affect the life of the implant
- Damaged instruments and devices must not be used for implant insertion
- The necks must not be used for trial articulation
- When using a cemented cup assembly, remove all loose bone cement; pieces of bone cement may penetrate between the friction surfaces and cause abnormal wear of the implant
- Always insert the implant into a precisely machined (verified by test components), clean and dry bed
- The articulation surfaces must be perfectly clean during the test and final insertion
- Clean and dry the anchoring and articulation surfaces thoroughly before inserting the definitive neck

After surgery

Patient

- The patient must be aware of the limitations caused by the implant
- Must be warned that his new implant can only be subjected to limited loading until the bone is fully healed
- In case of unexpected changes that may be related to the implant, it is recommended to see a specialist

Workplace performing surgery

- The load-bearing capacity of the implant cannot be compared to healthy bone
- Timely preventive checks of the artificial joint
- Based on the available information, MRI examinations with static magnetic field up to 3T can be performed in patients with implants from BEZNOSKA, s.r.o., but not earlier than 6 weeks after implantation and without signs of loosening.

Procedure during surgery

The trapeziometacarpal joint has its own detailed surgical procedure where the complete implantation is described.

Adverse complications

During surgery

- Damage to vascular and nervous structures
- Iatrogenic bone damage, sometimes leading to fracture
- Incorrect position of components, wrong resection

After surgery

- Cardiovascular disorders such as TEN, deep phlebotrombosis, MI, postoperative hematoma
- Wound healing disorders, infection
- Loosening of the entire total replacement or its component, deformity or fracture of the component

Product disposal

After use, the implant is listed as „hazardous waste“ from the healthcare sector. Dispose of according to the applicable legislation of the country where the implant is disposed of.

Manufacturer's recommendations:

The implant cannot be re-sterilised because of the risk:

- Damage to polished articulation surfaces, fasteners, anchoring surfaces
- Damage or contamination of the porous layer

Risks arising from repeated use from the patient's perspective:

- Failure of implants due to damage to articulating, connecting and anchoring surfaces or improper sterilization – risk of early reoperation
- Risk of infection in the patient – risk of difficult to treat infection and subsequent loosening
- Change in implant size or type, or incompleteness of the implant – risk of reoperation or early reoperation

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