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

www: www.beznoska.cz

EN Instructions for Use No. 001

Information on handling, possible applications and potential restrictions for the product group:

CEMENTED STEMS

Identification

Beznoška – Stem, 12/14, AK	Cemented revision stem, 12/14
Size range: 1/N – 3B	Size range: 11/160 – 14/300
CSC – Cemented stem, 12/14	TRIO – Cemented stem, 12/14
Size range: 0 – 5	Size range: 01 – 8
Femoral stem with partial femoral repl., 12/14	CSC – Centralizer
Size range: 1 – 5	Size range: 1 – 5
Trio – Centralizer	Femoral stem collar
Size range: 01 – 8	Size range: 4 – 10
Plug under cemented stem	Cerclage – Cerclage tape
Size range: D8 – D19	Size range: D8 – D19
Cerclage – Slider with pikes	Cerclage – Plate
	Size range: L170 – L230

Information about sterility

- The product is supplied by the manufacturer in a sterile pack
- Sterility of the implant in the protective sterile packaging is guaranteed for 5 years. Expiry dates are marked clearly on the protective sterile packaging
- Not to be used if the sterile packaging is damaged. The product cannot be resterilised

Intended use

Cemented acetabular stems are a part of a set intended for total hip replacement using bone cement. The stems are designed for implanting in a femur with the 12/14 neck intended for fitting a hip joint head. The remaining implants are auxiliary items.

Product description

Beznoška – Stem, 12/14, AK is supplied in versions using combinations of three different lengths and four different cross-sections (a total of 8 versions). The neck has a uniform length and cone – 12/14 (euro). The CCD angle is 140°. The stem is made of ISO 5832-9. The stem surface is polished, the neck cone has a structured surface.

Cemented revision stem, 12/14 is supplied in versions using combinations of four different diameters and eight stem lengths (32 versions). The neck has a uniform length and cone – 12/14 (euro). The CCD angle is 140°. The stem is made of ISO 5832-9. The stem surface is polished, the neck cone has a structured surface.

CSC – Cemented stem, 12/14 is supplied in six sizes. The neck has a uniform length and cone – 12/14 (euro). The CCD angle is 135°. The stem is made of ISO 5832-9. The stem surface is polished, the neck cone has a structured surface. The PMMA centralizer is connected to the lower part of the stem with a cone.

Trio – Cemented stem, 12/14 is supplied in ten sizes. The neck has a uniform length and cone – 12/14 (euro). The CCD angle is 135°. The stem is made of ISO 5832-9. The stem surface is polished, the neck cone has a structured surface. The PMMA centralizer is connected to the lower part of the stem with a cone.

Femoral stem with partial femoral repl., 12/14 is supplied in a version with a conical neck 12/14 (euro) for fitting a metal head. The stem is made of ISO 5832-9, the partial femur replacement and the distance ring are made of UHMWPE (ISO 5834-2), the collar of the 12/14 cone is made of ISO 5832-9. The stem surface is polished, the neck cone has a structured surface. Distance rings may be included with the stem with partial femur replacement, 12/14.

Femoral stem collars are made in four basic thicknesses of 4, 6, 8 and 10 mm. Each thickness has an opening sized 10, 11,5, 13 and 14,5 mm. The washer shape reflects the shape of the TEP stem collar (Beznoška – Stem, 12/14, AK) by BEZNOŠKA, s.r.o. The washers for the endoprosthesis collar are made of ultra-high molecular weight polyethylene UHMWPE (ISO 5834-2).

Plugs under cemented stem are supplied in six sizes with a diameter of 8, 10, 12, 14, 16 and 18 mm with an opening for the plug loader. The plugs under the stem come in a uniform length of 14 mm for all sizes. The plugs under the stem are made of ultra-high molecular weight polyethylene UHMWPE (ISO 5834-2).

Cerclage – the 7,5 mm wide and 0,5 mm thick cerclage tape has an eyelet at the end used for closing the loop and as a support for the tensioner while tightening the loop around the bone. The tape is glued into position around the bone with an opening for the loader on the other end. When placing the loop on structured parts of the bone, one or a maximum of two sliders with pikes are placed on the tape to guarantee the loop has a stable position. The set may be supplemented with a plate. A plate secured with the band provides greater stability of the assembly. The plate is supplied in lengths of 170, 200 and 230 mm. Its thickness is 2,5 mm and width 10 mm. The tape and the slider are made of pure titanium (ISO 5832-2), and the plate is made of wrought titanium alloy Ti6Al4V (ISO 5832-3). The surface finish is achieved by blasting.

Indications

- Adult population 18+, men, women
- Cases of primary and secondary coxarthrosis, idiopathic necrosis of the femoral head, inflammatory and post-inflammatory conditions (progressive polyarthritis and Bekhterev's disease), post-traumatic conditions (pseudoarthrosis of the femoral neck, acetabular fractures).

Contraindications

Carcinoma in the upper femur, internal contraindication, patient's inability to cooperate after the surgery, infections, documented allergy to or intolerance of metals (ISO 5832-9) and bone cement.

Combinations:

The cemented acetabular stems can be combined with all heads (Metal head, AK, Metal head, Co, Ceramic head (BILOX DELTA)) supplied by BEZNOŠKA, s.r.o. Trio – Centralizer – each of the sizes has its own centralizer. Femoral stem collar has the same opening as the Beznoška – Stem, 12/14, AK.

CSC – Centralizer:

CSC – Centralizer, PMMA	1	x	x				
	2			x			
	3				x		
	4					x	
	5						x

The products cannot be combined with products from other manufacturers.

Cautions and warnings

Before surgery

Patient

- The patient consents to the surgery and associated risks
- The patient must be mentally capable of understanding the meaning of the surgery and observing the postoperative regimen
- The patient must be aware of the restrictions caused by the implant

Workplace performing surgery

- Ensures highly sterile operating conditions
- Checks the sterility of the packaging (damage to the packaging and expiry dates)
- Tools for implantation must be complete and functional
- Ensures that damaged, non-sterile or reoperated implants are not used
- Patients receiving a hip joint replacement must be informed that the implant's lifetime depends on their weight and level of activity
- All necessary implantation components must be available
- Only ceramic heads supplied by BEZNOŠKA, s.r.o. may be used for sets with a ceramic head.

During surgery

- During application, the instruments supplied by the manufacturer must be used and the manufacturer's recommended surgical procedure must be followed
- Surgical gloves must be used when handling the product
- The precisely shaped conical neck part with a structured surface must be protected while handling the product. Any damage to this surface affects the precision of implanting and the lifetime of the head
- Testing heads must be used for testing the joint connection
- Always introduce the implant in a shaped bone marrow cavity
- The implant may only be combined with implants supplied by BEZNOŠKA, s.r.o. or implants expressly recommended for use by BEZNOŠKA, s.r.o.

After surgery

Patient

- The patient must be aware of the restrictions caused by the implant
- The patient must be warned that the new implant may only undergo limited stress until the bone has healed fully
- In the event of any unexpected changes that may relate to the implant, a visit to a specialist is recommended

Workplace performing surgery

- The load-bearing capacity of the implant is not comparable with that of a healthy bone
- Perform preventive checks of the artificial joint in a timely manner
- Based on the available information, magnetic resonance imaging with a static magnetic field up to 3T may be performed in patients with implants by BEZNOŠKA, s.r.o., if at least 6 weeks have passed from the date of implantation and if the implant does not display any signs of being loose.

Procedure during surgery

A specific detailed surgical procedure describing the entire implantation process applies to each of the cemented stems

- Manual – Beznoška – Stem, 12/14, AK
- Manual – Cemented revision stem
- Manual – CSC – Cemented stem
- Manual – Trio – Cemented stem
- Manual – Cerclage – Cerclage system

Adverse complications

During surgery

- Damage to vascular and nerve structures
- latrogenic damage to the bone, sometimes even the occurrence of fractures

After surgery

- Cardiovascular defects, such as TEN, deep phlebotrombosis, MI, postoperative hematoma
- Wound healing disorders, infection
- Loosening of the entire artificial joint or its components, luxation or subluxation of the artificial joint, deformation or breakage of a component, flexion contracture, limb shortening

Product disposal

After being used, the implant is classified as "dangerous waste" in health care. It is to be disposed of in accordance with the valid legislation of the country where the implant is being disposed of.

Risks arising from repeated use from the patient's perspective:

- Failure of the implant due to damage to the conical neck surface, anchoring surfaces (risk of rubbing, effect of rubbed parts and subsequent loosening of the implant) or inadequate sterilization, pain, potential need to reoperate
- Risk of infection in the patient (risk of infection and loosening of the implant, pain)
- Mistaken implant size or type, or incomplete – non-functioning implant, necessary to reoperate.

Manufacturer's recommendations:

Revision cemented stem (310905, 310915) – not to be used in patients with a body weight greater than 60 kg. CSC – Cemented stem (317000, 317001) – not to be used in patients with a body weight greater than 80 kg. TRIO – Cemented stem (321999, 322000, 322001) – not to be used in patients with a body weight greater than 80 kg. The implant cannot be resterilised. If the sterile packaging is damaged, do not use the medical device and send it back to the manufacturer. The manufacturer does not recommend combining the Metal head Co when revising with the original PE insert, which could be damaged in any way.

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