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
EN Instructions for Use No. 002

Information on handling, possible applications and potential restrictions for the product group:

Identification		
Metal head Co, 12/14		Metal head AK, 12/14
Size range: D22, D28, D32, D36		Size range: D28, D32, D36
Ceramic head (BIOLOX DELTA), 12/14		
Size range: D28, D32, D36		

Information about sterility

- The product is supplied by the manufacturer in a sterile pack
- Sterility of the implant in the protective sterile packaging is guaranteed for 5 years. Expiry dates are marked clearly on the protective sterile packaging
- Not to be used if the sterile packaging is damaged. The product cannot be resterilised

Intended use

Hip joint heads are part of a set intended for total hip replacement. The heads are designed for fitting on stems with a 12/14 cone.

Product description

The *Metal head Co* is supplied in a diameter of 22, 28, 32 and 36 mm; the 22 mm head is made with two different cone depths: L (long) and M (medium), the 28, 32 and 36 mm heads are made with five cone depths: XXL (extra long), XL (extended), L (long), M (medium) and S (short). The head is made of ISO 5832-12 material and the surface is polished to a mirror-like shine.

The *Metal head AK, 12/14* is supplied with a diameter of 28, 32 and 36 mm; have five cone depths: XXL (extra long), XL (extended), L (long), M (medium) and S (short). The head is made of ISO 5832-9 material and the surface is polished to a mirror-like shine.

The *Ceramic head (BIOLOX DELTA)* is supplied with a diameter of 28, 32 and 36 mm with three cone depths and neck length indication – L (long), M (medium) and S (short). The head is made of ISO 6474-2 material.

Indications

- Adult population 18+, men, women
- Cases of primary and secondary coxarthrosis, idiopathic necrosis of the femoral head, inflammatory and post-inflammatory conditions (progressive polyarthrits and Bekhterev's disease), post-traumatic conditions (pseudo-arthritis of the femoral neck, acetabular fractures...), certain bone tumours, etc.

Contraindications

Carcinoma in the upper femur, internal contraindication, patient's inability to cooperate after the surgery, infections, documented allergy to or intolerance of metals (ISO 5832-9, ISO 5832-12) and bone cement.

Combinations

The hip joint head can be combined with all stems (Beznoska – Stem, 12/14, AK, Cemented revision stem, CSC – Cemented system, 12/14, TRIO – Stem, SF – Stem, RMD – Cementless stem, Femoral stem with partial femoral repl., 12/14) supplied by BEZnosKA, s.r.o. **The products cannot be combined with products from other manufacturers.**

Caution and warnings

Before surgery

Patient

- The patient consents to the surgery and associated risks
- The patient must be mentally capable of understanding the meaning of the surgery and observing the postoperative regimen
- The patient must be aware of the restrictions caused by the implant

Workplace performing surgery

- Ensures highly sterile operating conditions
- Checks sterility of the packaging (damage to the packaging and expiry dates)
- Tools for implantation must be functional and complete
- Ensures that damaged, non-sterile or reoperated implants are not used
- Patients receiving a hip joint replacement must be informed that the implant's lifetime depends on their weight and level of activity
- The implant cannot be used in a set with a stem made of wrought stainless steel (ISO 5832-1).
- All necessary implantation components must be available

During surgery

- During application, the instruments supplied by the manufacturer must be used and the manufacturer's recommended surgical procedure must be followed
- Surgical gloves must be used when handling the product
- The polished spherical head surface must be protected, for example from falling or scratching, while handling the product. Any damage to this surface significantly reduces the lifetime of the articulated cartilaginous acetabulum surface
- Ensure perfect cleanliness of the spherical head surface and the articulated cartilaginous acetabulum surface.
- The precisely shaped conical neck part with a structured surface must be protected while handling the product. Any damage to this surface affects the precision implanting and the lifetime of the head.

After surgery

Patient

- The patient must be aware of the restrictions caused by the implant
- The patient must be warned that the new implant may only undergo limited stress until the bone has healed fully
- If there are any unexpected changes that may relate to the implant, a visit to a specialist is recommended

Workplace performing surgery

- The load-bearing capacity of the implant is not comparable with that of a healthy bone
- Perform preventive checks of the artificial joint in a timely manner
- Based on the available information, magnetic resonance imaging with a static magnetic field up to 3T may be performed in patients with implants by BEZnosKA, s.r.o., if at least 6 weeks have passed from the date of implantation and if the implant does not display any signs of being loose.

Procedure during surgery

Aspecific detailed surgical procedure describing the entire implantation process applies to each of the cemented stems.

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Adverse complications

During surgery

- Damage to vascular and nerve structures
- Iatrogenic damage to the bone, sometimes even occurrence of fractures

After surgery

- Cardiovascular defects, such as TEN, deep phlebotrombosis, MI, postoperative hematoma
- Wound healing disorders, infection
- Loosening of the entire artificial joint or its components, luxation or subluxation of the artificial joint, deformation or breakage of a component, flexion contracture, limb shortening

Product disposal

After being used, the implant is classified as “dangerous waste” in health care. It is to be disposed of in accordance with the valid legislation of the country where the implant is being disposed of.

Manufacturer's recommendations:

The implant cannot be resterilised. If the sterile packaging is damaged, do not use the medical device and send it back to the manufacturer. The manufacturer does not recommend combining the Metal head Co when revising with the original PE insert, which could be damaged in any way.

Risks arising from repeated use from the patient's perspective:

- Failure of the implant due to damage to the conical neck surface (imperfect connection between the stem and the head – risk of the head loosening leading to particles being rubbed between the stem cone and the head, metallosis caused by friction parts, potential need to reoperate), articulated surfaces (risk of rubbing, effect of rubbed parts and subsequent loosening of the implant) or inadequate sterilization, pain, potential need to reoperate
- Risk of infection in the patient (risk of infection, loosening of the implant or need to reoperate)
- Mistaken implant size or type, or incomplete – non-functioning implant, need to reoperate

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