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

EN Instructions for Use No. 003

Information on handling, possible applications and potential restrictions for the product group:

CEMENTED ACETABULAR CUPS

Identification

02 – Cemented cup	Beznoska – Cemented cup
Size range: D44 – D66	Size range: D44 – D64

Information about sterility

- The product is supplied by the manufacturer in a sterile pack
- Sterility of the implant in the protective sterile packaging is guaranteed for 5 years. Expiry dates are marked clearly on the protective sterile packaging
- Not to be used if the sterile packaging is damaged. The product cannot be resterilised

Intended use

Cemented cups are part of a set intended for total hip replacement using bone cement. The cups are intended for implanting in the acetabulum.

Product description

02 – Cemented cup is supplied in three standard versions – standard, with 10° offset and antiluxation type. The cup is manufactured in sizes D44, 46, 48, 50, 52, 54, 56, 58, 60 mm for D28 or 32 mm heads and in size range D54, 56, 58, 60, 62, 64, 66 mm for D36 mm heads. PMMA centralizers provide an even, two-millimetre layer of bone cement between the cup surface and the shaped acetabulum. The small perimeter collar helps to compress the cement while applying the cup. The cups are made of the cross-linked ultra-high molecular weight polyethylene UHMWPE (ISO 5834–2). A ring of wrought stainless steel ISO 5832–1 is fitted on the edge of the cup for easier X-ray identification.

Beznoska – Cemented cup is supplied in many sizes (for D28 and D32 mm heads) and versions (23 versions in total). The cups are made of ultra-high molecular weight polyethylene UHMWPE (ISO 5834–2). A ring of wrought stainless steel (ISO 5832–1) is fitted on the edge of the cup for easier X-ray identification.

Indications

- Adult population, 18+, men, women
- Cases of primary and secondary coxarthrosis, idiopathic necrosis of the femoral head, inflammatory and post-inflammatory conditions (progressive polyarthritis and Bekhterev’s disease), post-traumatic conditions (pseudo-arthritis of the femoral neck, acetabular fractures…) , certain bone tumours, etc.

Contraindications

Carcinoma in the upper femur, internal contraindication, patient’s inability to cooperate after the surgery, infections, documented allergy to or intolerance of metals (ISO 5832–1) and bone cement.

Combinations

The cemented acetabular cups can be combined with all heads (Metal head, AK, Metal head, Co, Ceramic head (BILOX DELTA)) supplied by the BEZNOSKA, s.r.o. company, provided the head diameter matches the cup’s inner diameter. **The products cannot be combined with products from other manufacturers.**

Caution and warnings

Before surgery

Patient

- The patient consents to the surgery and associated risks
- The patient must be mentally capable of understanding the meaning of the surgery and observing the postoperative regimen
- The patient must be aware of the restrictions caused by the implant

Workplace performing the surgery

- Ensures highly sterile operating conditions
- Checks sterility of the packaging (damage to the packaging and expiry dates)
- Tools for implantation must be functional and complete
- Ensures that damaged, non-sterile or reoperated implants are not used
- Patients receiving a hip joint replacement must be informed that the implant’s lifetime depends on their weight and level of activity
- The implant is not to be used in patients with a documented allergy to or intolerance of metals (ISO 5832–1), UHMWPE (ISO 5834–2) and bone cement
- All necessary implantation components must be available

During surgery

- During application, the instruments supplied by the manufacturer must be used and the manufacturer’s recommended surgical procedure must be followed
- Surgical gloves must be used when handling the product
- The precisely shaped articulated surface must be protected while handling the cup
- The implant must be always applied in a precisely shaped, clean and dry acetabulum evenly covered with bone cement
- The bone cement specifically intended for cemented cups must be used and the bone cement manufacturer’s instructions must be observed
- Ensure that the surface of tools used to introduce the implant or test the joint connection is not damaged. Any damage to the surfaces used to introduce the implant or to test the joint connection causes damage to the implant’s articulated surfaces and thus has a negative impact on the lifetime of the implant and the head
- Ensure all articulated surfaces are perfectly clean when testing or performing the final joint connection
- Remove any loose bone cement. Pieces of bone cement may spread between the contact surfaces and cause abnormal wear and tear on the implant

After surgery

Patient

- The patient must be aware of the restrictions caused by the implant
- The patient must be warned that the new implant may only undergo limited stress until the bone has fully healed
- If there are any unexpected changes that may relate to the implant, a visit to a specialist is recommended

Workplace performing surgery

- The load-bearing capacity of the implant is not comparable with that of a healthy bone
- Perform preventive checks of the artificial joint in a timely manner
- Based on the available information, magnetic resonance imaging with a static magnetic field up to 3T may

CE 1014

be performed in patients with implants by BEZNOSKA, s.r.o., if at least 6 weeks have passed from the date of implantation and if the implant does not display any signs of being loose.

Procedure during surgery

A specific detailed surgical procedure describing the entire implantation process applies to each of the cemented cups.

- Manuál – 02 – Cemented cup
- Manuál – Beznoska – Cemented cup

Adverse complications

During surgery

- Damage to vascular and nerve structures
- latrogenic damage to the bone, sometimes even the occurrence of fractures

After surgery

- Cardiovascular defects, such as TEN, deep phlebothrombosis, MI, postoperative hematoma,
- Wound healing disorders, infection
- Loosening of the entire artificial joint or its components, luxation or subluxation of the artificial joint, deformation or breakage of a component, flexion contracture, limb shortening

Product disposal

After being used, the implant is classified as “dangerous waste” in health care. It is to be disposed of in accordance with the valid legislation of the country where the implant is being disposed of.

Risks arising from repeated use from the patient’s perspective:

- Failure of the implant due to damage to the articulated surface (loosening due to friction on the parts), anchoring surfaces (risk of rubbing, effect of rubbed parts and subsequent loosening of the implant) or inadequate sterilization, pain, potential need to reoperate
- Risk of infection in the patient (risk of infection and loosening of the implant, pain or need to reoperate)
- Mistaken implant size or type, or incomplete – non-functioning implant, necessary to reoperate

Manufacturer’s recommendations:

The implant cannot be resterilised. If the sterile packaging is damaged, do not use the medical device and send it back to the manufacturer. The manufacturer does not recommend combining the Metal head Co when revising with the original PE insert, which could be damaged in any way.

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