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CE **1014**

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EN **Instructions for Use No. 010**

CE **1014**

Information on handling, possible applications

and potential restrictions for the product group:

CERVICAL-CAPITAL ENDOPROSTHESIS

Beznoska – Stem CEP	Metal head, 12/14
Size range: D38 – D56	Size range: D38 – D56
Bi-articular head, 12/14	CSB – Stem
Size range: D42 – D56	Size range: D40 – D56
CSB – Stem, 12/14	
Centralizer, 01 – 8	CSC – Centralizer, 2

Information about sterility

- The product is supplied by the manufacturer in a sterile pack
- Sterility of the implant in the protective sterile packaging is guaranteed for 5 years. Expiry dates are marked clearly on the protective sterile packaging
- Not to be used if the sterile packaging is damaged. The product cannot be resterilised

Intended use

The implants are designed as a cervical-capital hip joint replacement using bone cement.

Product description

Beznoska – Stem CEP is intended for implanting in the femur using bone cement. It is supplied with a uniform stem and various diameters of the hollow head in a total of 10 versions from D38 to 56 mm. CEP stems are made of stainless steel (ISO 5832-1). The endoprosthesis has a polished surface and the head surface is polished to a mirror-like shine.

Metal head, 12/14, is intended for fitting on cemented stems with a conical neck 12/14. It is produced in 10 sizes from D38 to 56 mm. The Metal head, 12/14, is made of wrought stainless steel (ISO 5832-1). The head surface is polished to a mirror-like shine.

Bi-articular head, 12/14, is intended for fitting on cemented TEP stems with a conical neck 12/14. It is produced in sizes from D42 to 56 mm. The implant comes as a fixed set and cannot be dismantled. The head cup is made of wrought stainless steel (ISO 5832-1), the inner metal head is made of wrought high-nitrogen stainless steel (ISO 5832-9). The outer surfaces of both parts are polished and the inner head surface is polished to a mirror-like shine. The lining of the head cup is made of ultra-high molecular weight polyethylene UHMWPE (ISO 5834-2).

CSB – Stem is intended for implanting in the femur using a PMMA centralizer. CSB – Stem is supplied in nine head sizes from D40 to 56 mm. CSB – Stem, 12/14, is supplied in a single version and the endoprosthesis neck is fitted with a collar. The CCD angle is 135°. The stems are made of wrought high-nitrogen stainless steel (ISO 5832-9), the head is made of wrought stainless steel (ISO 5832-1). The stem surface is polished. The PMMA centralizer is connected to the lower stem part with a cone.

Indications

- Adult population 18+, men, women
- Cases are indicated in primary and secondary coxarthrosis, idiopathic head necroses, fractures and post-traumatic conditions of the femoral neck in patients of older biological age.

Contraindications

Biologically younger or biologically unsuitable patients, where full joint replacement can be considered, significantly dysplastic acetabulum, osteoporosis, internal contraindications, infection, major osteoarthritis.

Combinations

The cervical-capital stems with a 12/14 cone can be combined with all heads (Metal head, AK, Metal head, Co, Ceramic head (BILOX DELTA)) supplied by BEZŇOSKA, s.r.o. The centralizer is suitable for all CSB stems. **The products cannot be combined with products from other manufacturers.**

Caution and warnings

Before surgery

Patient

- The patient consents to the surgery and associated risks
- The patient must be mentally capable of understanding the significance of prevention, as well as observing prevention.
- The patient must be aware of the restrictions caused by the implant

Workplace performing surgery

- Ensure highly sterile operating conditions
- Check the sterile packaging for potential damage and the expiry dates on the packaging
- Tools for implantation must be complete and functional
- The implant is intended for single use only
- Damage or reoperated implants must not be used under any circumstances
- The implant is not to be used in patients with a documented allergy to or intolerance of metals (ISO 5832-1 alloy) and bone cement.

During surgery

- During application, the instruments supplied by the manufacturer must be used and the manufacturer's recommended surgical procedure must be followed
- The polished spherical head surface must be protected, for example from falling or scratching, while handling the product. Any damage to this surface significantly reduces the lifetime of the articulated cartilaginous acetabulum surface.
- Ensure perfect cleanliness of the spherical head surface and the articulated cartilaginous acetabulum surface.
- Always introduce the implant in a shaped bone marrow cavity.

After surgery

Patient

- The patient must be aware of the restrictions caused by the implant
- The patient must be warned that the new implant may only undergo limited stress until the bone has healed fully
- In the event of any unexpected changes that may relate to the implant, a visit to a specialist is recommended

Workplace performing surgery

- The load-bearing capacity of the implant is not comparable with that of a healthy bone
- Perform preventive checks of the artificial joint in a timely manner
- Based on the available information, magnetic resonance imaging with a static magnetic field up to 3T may be performed in patients with implants by BEZŇOSKA, s.r.o., if at least 6 weeks have passed from the date of implantation and if the implant does not display any signs of being loose.

Procedure during surgery

A specific detailed surgical procedure describing the entire implantation process applies to each of the cervical-capital endoprotheses.

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Adverse complications

During surgery

- Damage to vascular and nerve structures
- latrogenic damage to the bone, sometimes even occurrence of fractures

After surgery

- Cardiovascular defects, such as TEN, deep phlebothrombosis, MI, postoperative hematoma
- Wound healing disorders, infection
- Loosening of the entire CEP, luxation or subluxation of the CEP, deformation or breakage of a component, flexion contracture, limb shortening

Product disposal

After being used, the implant is classified as “dangerous waste” in health care. It is to be disposed of in accordance with the valid legislation of the country where the implant is being disposed of.

Manufacturer's recommendations:

The implant cannot be resterilised due to the risk of:

- Damage to the anchoring surfaces
- Damage to the polished articulated surface

Risks arising from repeated use from the patient's perspective:

- Failure of the implant due to damage to the anchoring surfaces or inadequate sterilization, pain, potential need to reoperate
- Risk of patient developing infection
- Mistaken implant size or type – non-functioning implant, necessary to reoperate

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