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
CE 1014

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
EN Instructions for Use No. 013

CE 1014

Information on handling, possible applications
and potential restrictions for the product group:

AUXILIARY ITEMS

Identification

K-wire	Osteofix wire
Size range: D1.2 – D3, L100-L300	Size range: D0.5 – D1.2

Information about sterility

- The product is supplied by the manufacturer in a non-sterile condition
- The sterilisation method is specified in the attached sterilisation procedure for non-sterile implants no. 500
- The product can be resterilised

Intended use

The K-wire or Osteofix wire are auxiliary items used for a total hip replacement.

Product description

KWire is a product with a round cross-section, a tip (triangular or pike) with cutting edges to make introduction easier and an end piece for fixing into a head or a chuck. The product is characterised by the diameter of the cylindrical section and the total length. K-wires are supplied in sizes with a diameter from 1.2 to 3 mm and a length of 100 to 300 mm. The wires are made of wrought stainless steel (ISO 5832-1). The product can also be used as an implant. The surface finish is achieved by polishing.

Osteofix wire is used for cerclage, fixing plates supplied by BEZNOSKA, s.r.o. for pectus excavatum, for auxiliary fixation or as a diagnostic wire for transparent implants. Osteofix wire is supplied in rolls containing 5 m of the wire with diameters of 0.5, 0.6, 0.7, 0.8, 1.0 and 1.2 mm. It is made of wrought stainless steel (ISO 5832-1). The surface finish is achieved by polishing.

Indications

- Adult population, 18+, men, women
- Indications are auxiliary fixation, temporary or main fixation of bone fragments.

Contraindications

Internal contraindication, infections, documented allergy to or intolerance of metals (ISO 5832-1). The product does not meet the requirements for fully stable osteosynthesis.

Combinations:

The implants can be combined exclusively with implants from BEZNOSKA, s.r.o. suitable for the relevant purpose.

The products cannot be combined with products from other manufacturers.

Caution and warnings

Before surgery

Patient

- The patient consents to the surgery and associated risks
- The patient must be mentally capable of understanding the meaning of the surgery and observing the postoperative regimen
- The patient must be aware of the restrictions caused by the implant

Workplace performing surgery

- Ensures highly sterile operating conditions
- Tools for implantation must be functional and complete
- The implant is only intended for single use
- Ensures that damaged, non-sterile or reoperated implants are not used
- The implant is not to be used in patients with a documented allergy to or intolerance of metals (alloy ISO 5832-1)

During surgery

- During application, the instruments supplied by the manufacturer must be used and the manufacturer's recommended surgical procedure must be followed
- Surgical gloves must be used when handling the product
- The wires can be shaped as needed using regular tools
- While handling the product, caution must be taken to prevent any damage to the working parts (sharp tips and ends) and bending the cylindrical part of the wire, when used as a guide wire

After surgery

Patient

- The patient must be aware of the restrictions caused by the implant
- The patient must be warned that the new implant may only undergo limited stress until the bone has fully healed
- If there are any unexpected changes that may relate to the implant, a visit to a specialist is recommended

Workplace performing surgery

- The load-bearing capacity of the implant is not comparable with that of a healthy bone
- Perform preventive checks of the artificial joint in a timely manner
- Based on the available information, magnetic resonance imaging with a static magnetic field up to 3T may be performed in patients with implants by BEZNOSKA, s.r.o., if at least 6 weeks have passed from the date of implantation and if the implant does not display any signs of being loose.

Procedure during surgery

These auxiliary items are intended for use by experienced surgeons and their use is entirely at the surgeon's discretion

Adverse complications

During surgery

- Damage to vascular and nerve structures
- Iatrogenic damage to the bone, sometimes even the occurrence of fractures

After surgery

- Cardiovascular defects, such as TEN, deep phlebothrombosis, MI, postoperative hematoma
- Wound healing disorders, infection
- Loosening of the entire artificial joint or its components, luxation or subluxation of the artificial joint, deformation or breakage of a component, flexion contracture, limb shortening

Product disposal

After being used, the implant is classified as "dangerous waste" in health care. It is to be disposed of in accordance with the valid legislation of the country where the implant is being disposed of.

Manufacturer's recommendations:

The implant can be resterilised:

- Observe the sterilisation instructions provided by the manufacturer

Risks arising from repeated use from the patient's perspective:

- Failure of the implants due to damage to the articulated surface, anchoring surfaces due to damage to the anchoring and fixation elements – risk of rubbing, breakage of fixation screws or inappropriate sterilization, pain or need to reoperate
- Risk of infection in the patient (risk of infection and loosening of the implant, pain)
- Mistaken implant size or type, or incomplete – non-functioning implant, need to reoperate

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