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

Information on handling, possible applications and potential restrictions for the product group:

## TOTAL KNEE REPLACEMENT – SVS

### Identification

SVS – Femoral component	SVL – Tibial plateau
Size range: 1 – 6	Size range: 1 – 6
SVS – Insert	Patella
Size range: 1 – 6	Size range: D25 – D34

### Information about sterility


- The product is supplied by the manufacturer in a sterile pack
- Sterility of the implant in the protective sterile packaging is guaranteed for 5 years. Expiry dates are marked clearly on the protective sterile packaging
- Not to be used if the sterile packaging is damaged. The product cannot be re-sterilised

### Intended use

The SVS cemented implant is designed as a total knee replacement using bone cement.

### Product description

The stabilized cemented knee replacement is designed for primary knee joint replacement using bone cement during surgeries with excision of the anterior cruciate ligament. The replacement set comprises of the SVS type femoral component the SVL type tibial plateau, the SVS type insert and the connecting screw. It may be supplemented with the patella replacement. The femoral and tibial parts and patella sizes (but not versions) may be combined (see the table).

		SVL/N – Tibial plateau and SVS – Insert					Patella				
		1	2	3	4	5	6	1	2	3	4
SVS – Femoral component	2										
	3		2					X			
	4		3	X				X	X	X	
	4		4/3		X			X	X	X	X
	5		4			X		X	X	X	X
	6		5				X	X	X	X	X
		6					X	X	X	X	

The table always applies to the R (right) or L (left) versions

### SVS – Femoral component

The stabilized femoral component is designed for total knee joint replacement with excision of the anterior cruciate ligament and is intended for application using bone cement. It is made in 6 sizes in L/R versions. The articulated part of the component comprises of two condyles joined by a front disc. It is important to create a limited intercondylar space for guiding the elevated nose of the insert allowing excision of the cruciate ligament. The component is made of cobalt alloy (ISO 5832-4). The articulated surfaces are polished, and the internal surfaces are finely blasted.

### SVL – Tibial plateau

The standard tibial plateau is designed for total knee joint replacement and for application using bone cement. It is made in 6 sizes in left and right versions. The component includes an upper plate designed for fitting on a precisely resected proximal part of the tibia and the stem with ribs. The plate has an oval, asymmetric shape (left and right) and has an opening for the anterior cruciate ligament. The plate has an oval depression for the insert on the upper side and an opening for the screw. The opening is designed for the connecting screw securing the insert. The lower part of the plate is structured with recesses designed for primary fixation of the component in bone cement. The component is anchored in the tibia in particular with a conical stem and a pair of anti-rotational ribs. The component is made of wrought titanium alloy Ti6Al4V (ISO 5832-3), which has been used by leading global manufacturers for a long time. The surface of the implant is finely blasted.

### SVS – Insert

The stabilized insert is designed for total knee joint replacement with excision of the anterior cruciate ligament. It is made in 6 sizes and 5 thicknesses (8, 10, 12, 15 and 18 mm) and in left and right versions. The insert has an oval, asymmetric shape and has a cut-out and a lighter section. The implant's articulated surface comprises of two cylindrical surfaces divided by a bridge to ensure the replacement functions following excision of the anterior cruciate ligament, the femoral component has frontal stability and there is limited rotation around the vertical axis. The bridge has an opening for the connecting screw to ensure a secure connection with the tibial plateau (the screw is a standard part of the insert pack). The insert is made of ultra-high molecular weight polyethylene (ISO 5834-2) and the screw is made of wrought titanium alloy (ISO 5832-3). The surface of the insert is polished, and the connecting screw surface is finely blasted.

### Patella

The patella replacement is designed for total knee joint replacement with the option of maintaining the anterior cruciate ligament and for replacement with anterior stabilization. It is made in four sizes. The component is rotationally symmetric with a concave articulated surface. The anchoring part comprises of a cylindrical fitting with crosswise grooves preventing rotation in the mounting. All sizes are made with identical anchoring. The component is made of ultra-high molecular weight polyethylene (ISO 5834-2). The surface of the implant on the articulation surface is very finely machined.

### Indications

- Adult population 18+, men, women
- Severe arthropathy based on rheumatoid arthritis or osteoarthritis, severe post-traumatic conditions and certain systemic diseases. The patella replacement is indicated for major degenerative changes of the patellar cartilage and patella tracking disorders. The patella replacement should be performed preferentially in the event of TEP primary implantation.

### Contraindications

Infectious diseases or local infections, severe neuromuscular or vascular diseases. Poor quality of bone structures, severe knee instability in the area of collateral ligaments. The implant is not to be used in patients with a documented allergy to or intolerance of metals (ISO 5832-3, ISO 5832-4), UHMWPE (ISO 5834-2) and bone cement. The patella replacement is contraindicated in cases of patellar osteoporosis and small patella thickness and size.

### Caution and warnings

#### Before surgery

##### Patient

- The patient consents to the surgery and associated risks
- The patient must be mentally capable of understanding the meaning of the surgery and observing the postoperative regimen
- The patient must be aware of the restrictions caused by the implant

#### Workplace performing surgery

- Ensures highly sterile operating conditions
- Checks sterility of the packaging (damage to the packaging and expiry dates)
- Tools for implantation must be complete and functional
- The implant is intended for single use only
- Ensures that damaged, non-sterile or reoperated implants are not used
- Patients receiving a knee joint replacement must be informed that the implant's lifetime depends on their weight and level of activity
- All necessary implantation components must be available

#### During surgery

- During application, the instruments supplied by the manufacturer must be used and the manufacturer's recommended surgical procedure provided in the manual must be followed
- Surgical gloves must be used when handling the product
- The polished articulated surfaces, the surfaces designed for fitting the insert and all threaded surfaces must be protected when handling the product. Any damage to the articulated surfaces adversely affects the lifetime of the implant
- Damaged tools must not be used to introduce the insert
- The SVS – Insert must not be used for testing the joint connection
- Do not use the SVS – Femoral component to test the joint connection
- The SVL – Tibial plateau can be used for testing the joint connection after cementing (in a set with a test articulation insert)
- Remove any loose bone cement; pieces of bone cement may spread between the contact surfaces and cause abnormal wear and tear on the implant
- Ensure all articulated surfaces are perfectly clean when testing or performing the final joint connection
- The Patella must not be used for testing the joint connection
- The upper surface of the SVL – Tibial plateau must be perfectly clean and dry before implanting the SVS – Insert – rinse and dry the surfaces thoroughly before placing the insert
- The connecting screw must be tightened appropriately (2 to 2.5 Nm)
- The implant may only be combined with implants supplied by BEZNOSKA, s.r.o. or implants expressly recommended for use by BEZNOSKA, s.r.o.

#### After surgery

##### Patient

- The patient must be aware of the restrictions caused by the implant
- The patient must be warned that the new implant may only undergo limited stress until the bone has healed fully
- In the event of any unexpected changes that may relate to the implant, a visit to a specialist is recommended

#### Workplace performing surgery

- The load-bearing capacity of the implant is not comparable with that of a healthy bone
- Perform preventive checks of the artificial joint in a timely manner
- Based on the available information, magnetic resonance imaging with a static magnetic field up to 3T may be performed in patients with implants by BEZNOSKA, s.r.o., if at least 6 weeks have passed from the date of implantation and if the implant does not display any signs of being loose.

#### Procedure during surgery

A detailed surgical procedure describing the complete implantation process is available for the total knee replacement – SVS type.

#### Adverse complications

##### During surgery

- Damage to vascular and nerve structures
- Iatrogenic damage to the bone, sometimes even fractures may occur
- Incorrect position of the components
- Incorrect resection

##### After surgery

- Cardiovascular defects, such as TEN, deep phlebotrombosis, MI, postoperative hematoma
- Wound healing disorders, infection
- Loosening of the entire total replacement or its components, deformation or breakage of a component, flexion contracture, limb shortening

#### Product disposal

After being used, the implant is classified as "dangerous waste" in health care. It is to be disposed of in accordance with the valid legislation of the country where the implant is being disposed of.

#### Manufacturer's recommendations:

The implant cannot be re-sterilised due to the risk of:

- Tibial component – damage to the conical and anchoring surfaces (rubbing of the cement and stem)
- Insert – damage to the articulated surface (release of rubbed parts), damage to insertion elements, material degradation (high risk of component destruction)
- Femoral component – damage to the articulated surface (release of rubbed parts), damage to the anchoring elements

#### Risks arising from repeated use from the patient's perspective:

- Failure of the implants due to damage to the articulated and conical surface, anchoring surfaces (risk of rubbing, effect of rubbed parts and subsequent loosening of the implant), pain or inadequate sterilization – risk of premature need to reoperate
- Risk of infection in the patient (risk of infection and loosening of the implant) – risk of urgent need to reoperate
- Mistaken implant size or type or incomplete implant – risk of repeated surgery or premature need to reoperate

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