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

Information on handling, possible applications and potential restrictions for the product group:

## TOTAL KNEE REPLACEMENT – SVR

### Identification

SVR – Femoral component	SVR – Femoral coupler
SVR – Cementless femoral stem	SVR – Distal femoral aug.
SVR – Dorsal femoral aug.	SVR II – Fixed tibial component
SVR II – Rotation tibial component	SVR II – Offset
SVR II – Fixed insert	SVR II – Rotary insert
SVR – Cementless tibial stem	SVR – Shaped tibial aug.
SVR – Block tibial aug.	SVR II – Metaphyseal tibial aug.
SVR II – Stem for metaphyseal tibial aug.	SVR II – Stem plug for metaphyseal tibial aug.
SVR II – Femoral component for metaphyseal aug.	SVR II – Metaphyseal femoral aug.
SVR II – Stem for metaphyseal femoral aug.	SVR I – Coupling stem for femoral metaphyseal aug.
SVR II – Sleeve for metaphyseal coupling stem	Patella

### Information about sterility


- The product is supplied by the manufacturer in a sterile pack
- Sterility of the implant in the protective sterile packaging is guaranteed for 5 years. Expiry dates are marked clearly on the protective sterile packaging
- Not to be used if the sterile packaging is damaged. The product cannot be re-sterilised

### Intended use

The implants are part of a total knee replacement set with or without bone cement.

### Product description

The revision knee replacement is designed for primary and revision knee joint replacement using bone cement during surgeries with excision of the anterior cruciate ligament. The replacement set always comprises of a type SVR or SVR II femoral component, a type SVR II tibial component, a type SVR insert and extended stems. The set may be supplemented with augmentations or a patella replacement as needed. The femoral and tibial parts and patella sizes may be combined to a limited extent (see the table).

### SVR – Femoral component

#### SVR II – Femoral component for metaphyseal aug.

The SVR and SVR II femoral components are designed for total knee joint replacement with excision of the anterior cruciate ligament and are intended for application using bone cement. They are made in five sizes in left and right versions. The articulation surfaces of the condyles and the external surfaces are based on the anatomic conditions of the knee joint to ensure the total replacement functions optimally. The intercondylar space for guiding the elevated nose of the insert is an important element allowing excision of the cruciate ligament and the substitution of its function. The internal surfaces are structured with recesses to securely fit the component into

bone cement. The precisely manufactured profile of the internal surfaces and the intercondylar space ensures the component fits precisely and has a secure connection with the bone. The SVR and SVR II types differ in a major structural component – connection with an extended femoral stem. The SVR type has grooves and a pair of spherical recesses on the upper side for fitting the coupler and the extended stem and setting in two positions in the A-P direction. The SVR II type, intended exclusively for combination with metaphyseal augmentation, has a fixed conical extension on the upper side of the structure for fitting the stem and metaphyseal augmentation set. The extension does not allow movement in the A-P direction or selection of the Q-angle, which is fixed (5°). The cylindrical countersinking on the inner sides of the condyles is designed for augmentation hinges. The component is made of cobalt alloy (ISO 5832-4). The articulated surfaces are polished, and the internal surfaces are finely blasted.

#### SVR II – Stem for metaphyseal femoral aug.

The stem is used as an auxiliary guiding element for fitting the femoral component set in the marrow cavity and for a precise and solid connection of the femoral component with the metaphyseal augmentation. The stem is connected to the femoral component through a cone and a screw. The outer conical surface is used for fitting the set in the metaphyseal augmentation. The stems are made in 3 diameters (12, 14 and 16 mm) and lengths 150 and 210 mm. The range of stems includes a transition stem designed for use in combination with the sleeve to provide additional cement fixation of the set in the marrow channel of the femur. The transition stem has a uniform diameter of 12 mm and two lengths (150 and 210 mm). The stem is made of titanium alloy (ISO 5832-3). The surface of the implant is finely blasted.

#### SVR II – Sleeve for metaphyseal coupling stem

The sleeve is used for additionally fixing the femoral component set stem into the marrow cavity. The sleeve is connected to the transition stem through a cone. The external surface of the cover is grooved similarly to the standard tibial extension stems. The covers are made in five diameters (14, 16, 18, 20 and 22 mm) and a uniform length of 55 mm. The sleeve is made of titanium alloy (ISO 5832-3). The surface of the implant is finely blasted.

#### SVR II – Metaphyseal femoral aug.

The augmentation is designed for assembling with the SVR II type femoral component. It is intended for application without bone cement. The augmentation is made in five sizes. The augmentation has a stepped outer shape derived from the shape of the cavity of the metaphyseal part of the femur. The two upper steps are treated with a Ti spray coat (ISO 5832-2) and the remaining external surfaces are modified with nanostructuring. The internal conical opening is intended for connecting to the femoral stem for metaphyseal augmentations – the joint is self-locking. The augmentation is made of titanium alloy (ISO 5832-3).

#### SVR II – Fixed tibial component

The component is designed for application using bone cement. It is made in six sizes (1–6). The component includes an upper plate, designed for fitting on a precisely resected proximal part of the tibia, and a stem with ribs. The plate has a symmetric oval shape and a cut-out in the rear part. The plate has an oval depression for the insert on the upper side and a conical opening for the fixing stem. The lower part of the plate is structured with recesses designed to fix the component in bone cement, for example for inserting the augmentations. Two pairs of openings with an M6 thread for securing augmentations with screws are found in the recesses – the openings are sealed with PE plugs in the delivered component. The component has a pair of anti-rotational ribs in the transition between the plate and the stem. The conical component stem, in combination with other elements, ensures the stability of the implant. It has a conical surface for connecting the offset or the extended stem of metaphyseal augmentations onto the outer surface of the lower part and an opening for the extended stem securing screw. The component is made of titanium alloy (ISO 5832-3). The surface of the implant is finely blasted.

#### SVR II – Rotation tibial component

The component is designed for application using bone cement. It is made in five sizes (2–6). The component includes an upper plate, designed for fitting on a precisely resected proximal part of the tibia, and a stem with ribs. The plate has a symmetrical oval shape and a cut-out in the rear part. The upper surface of the plate is perfectly flat and polished for articulation with the rotational insert and is fitted with a conical opening for the insert fixation stem. The polished surface has protrusions on the dorsal side for limiting insert rotation. As with the fixed tibial component, the bottom part of the plate is structured with grooves for fixing the component into bone cement or for inserting augmentations and includes openings with an M6 thread for augmentation screws. The component also has a pair of anti-rotational ribs. The conical component stem, in combination with other elements, ensures the stability of the implant. It has a conical surface for connecting the offset or the extended stem of metaphyseal augmentations on the outer surface of the lower part and an opening for the extended stem securing screw. The component is made of cobalt alloy (ISO 5832-4). The surface of the implant is polished or finely blasted.

#### SVR II – Fixed insert

The insert is designed for a set with a fixed type of the tibial component. It is made in six sizes and six thicknesses (10, 12, 15, 18, 22 and 26 mm). The insert has a symmetric oval shape. The articulation surface of the implant comprises of two cylindrical surfaces divided by a bridge (an elevated "intercondylar eminence") to ensure the replacement functions properly following excision of the anterior cruciate ligament, and to ensure frontal stability of the femoral component and limited rotation around the vertical axis. The insert has a fixing stem in its lower part. A solid connection with the tibial component is provided by a securing screw. The insert is made of UHMWPE (ISO 5834-2). The insert surface is finely machined.

#### SVR II – Rotary insert

The insert is designed for a set with the rotational type of tibial component. It is made in five sizes and six thicknesses (10, 12, 15, 18, 22 and 26 mm). The insert has a symmetric oval shape. The articulation part of the insert comprises of a shaped contact surface for the femoral component. The surfaces are divided by a bridge (elevated "intercondylar eminence") to ensure the replacement functions properly following excision of the anterior cruciate ligament, and to ensure frontal stability of the femoral component and limited rotation around the vertical axis. A conical stem ensuring a rotational connection with the tibial component and a pair of recesses limiting rotation to the extent of  $\pm 10^\circ$  are in the lower part of the insert. The insert is made of UHMWPE (ISO 5834-2). The insert surface is finely machined.

#### SVR II – Offset

The offset is used as a transition piece between the tibial component and the extended stem allowing the extended stem axis to be shifted in relation to the axis of the conical stem of the tibial component. The shift of the axis can be set in planes passing through the tibial stem axis in 30° turns; the offset can be tuned in relation to the axis by 0, 3 and 6 mm. The offset is connected to the tibial component via a conical joint secured with a screw. The extended stem is screwed into the offset and is inserted into the offset in a self-locking cone. The offset is made of titanium alloy (ISO 5832-3). The surface of the implant is finely blasted.

#### SVR II – Stem for metaphyseal tibial aug.

The stem is used as an auxiliary guiding element for fitting the tibial component set in the marrow cavity and for a precise and solid connection of the tibial component with the metaphyseal augmentation. The stem is connected to the tibial component in the same manner as the offset (a cone with a screw). The outer conical surface is used for fitting the set into the metaphyseal augmentation. The stems are made in four diameters (10, 12, 14 and 16 mm) and lengths of 100, 120 and 150 mm. The range of stems includes a short stem plug for sealing the opening in the metaphyseal augmentation if the use of a standard-length stem is not appropriate or possible. The stem is made of titanium alloy (ISO 5832-3). The surface of the implant is finely blasted.

#### SVR II – Metaphyseal tibial aug.

The augmentation is designed mainly for assembly with the tibial component in the rotational version. It is intended for application without bone cement. The augmentation is made in four sizes and two versions (for the tibial component 1-3 version A and 4-6 version B). The augmentation has a stepped outer shape derived from the

shape of the cavity for the metaphyseal part of the tibia. The two upper steps are treated with a Ti spray coat and the remaining external surfaces are modified with nanostructuring. The internal conical opening is intended for connecting to the tibial stem for metaphyseal augmentations – the joint is self-locking. The shaped recesses in the inner part allow the augmentation to turn in relation to the tibial component by  $\pm 15^\circ$  (1A, 1B types),  $\pm 20^\circ$  (2B – 4B, 2A – 4A types). The augmentation is made of titanium alloy (ISO 5832-3). The implant surface is covered with a nanostructure and with a spray coat (ISO 5832-2) in the upper part.

#### Patella

The patella replacement is designed for total knee joint replacement with the option of maintaining the anterior cruciate ligament and for replacement with anterior stabilization. It is made in four sizes identical for the left and right versions. The component is rotationally symmetric with a convex articulation surface. The anchoring part comprises of a cylindrical fitting with crosswise grooves preventing rotation in the mounting. All sizes are made with identical anchoring. The component is made of UHMWPE (ISO 5834-2). The surface of the implant on the articulation surface is very finely machined.

#### Related implants that can be used in a revision replacement set

##### SVR – Femoral coupler

SVR – Cementless femoral stem

SVR – Cementless tibial stem

SVR – Femoral augmentation (distal, dorsal)

SVR – Block augmentation (shaped, block)

### Indications

- Adult population 18+, men, women
- Failures caused by ligament instability or extensor deficiency, aseptic loosening of the components, conditionally also complications in the patellofemoral joint, osteolysis caused by polyethylene rubbing, post-traumatic cases (periprosthetic fractures).

### Contraindications

Infectious diseases or local infections, severe neuromuscular or vascular diseases, poor quality bone structures, severe knee instability in the area of collateral ligaments, documented allergy to or intolerance of metals (ISO 5832-4, ISO 5832-3, ISO 5832-2), UHMWPE (ISO 5834-2) and bone cement. The patellar replacement is contraindicated in osteoporosis of the patella and in the event of small patella thickness or size.

### Caution and warnings

#### Before surgery

##### Patient

- The patient consents to the surgery and associated risks
- The patient must be mentally capable of understanding the significance of prevention, as well as observing prevention.
- The patient must be aware of the restrictions caused by the implant

##### Workplace performing surgery

- Ensures highly sterile operating conditions
- Checks the sterility of the packaging (damage to the packaging and expiry dates)
- Tools for implantation must be complete and functional
- Ensures that damaged, non-sterile or reoperated implants are not used
- Ensures the completeness of the implant set, instruments and equipment required for implantation in exceptional cases, decides on its use in cases of severe osteoporosis, suffered infections, significantly overweight patients, patients with major physical stress and patients addicted to drugs or alcohol
- Patients receiving a knee joint replacement must be informed that the implant's lifetime depends on their weight and level of activity

### During surgery

- During application, the instruments supplied by the manufacturer must be used and the manufacturer's recommended surgical procedure provided in the printed manual or the instructional video must be followed
- The polished articulation surfaces, the surfaces designed for fitting the insert and all threaded surfaces must be protected when handling the product. Any damage to the articulation surfaces adversely affects the lifetime of the implant
- Damaged tools must not be used to introduce the insert
- The insert must not be used for testing the joint connection
- The femoral component and the tibial component may only be used to test the joint connection after cementing (in the set with a test articulation insert)
- Remove any loose bone cement; pieces of bone cement may spread between the contact surfaces and cause abnormal wear and tear on the implant
- Ensure all articulated surfaces are perfectly clean when testing or performing the final joint connection
- The patella replacement must not be used for testing the joint connection
- The upper surface of the tibial component must be perfectly clean and dry before implanting the insert – rinse and dry the surfaces thoroughly before placing the insert
- The connecting screws and other threaded joints must be tightened appropriately

### After surgery

##### Patient

- The patient must be aware of the restrictions caused by the implant
- The patient must be warned that the new implant may only endure limited stress and that its load-bearing capacity increases gradually and individually – according to the physician's instructions
- In the event of any unexpected changes that may relate to the implant, a visit to a specialist is recommended

##### Workplace performing surgery

- The load-bearing capacity of the implant is not comparable with that of a healthy bone
- Perform preventive checks of the artificial joint in a timely manner
- Based on the available information, magnetic resonance imaging with a static magnetic field up to 3T may be performed in patients with implants by BEZnoska, s.r.o., if at least 6 weeks have passed from the date of implantation and if the implant does not display any signs of being loose.

### Procedure during surgery

A detailed surgical procedure describing the complete implantation process is available for the total knee replacement – SVR type.

### Adverse complications

#### During surgery

- Damage to vascular and nerve structures
- Iatrogenic: damage to the bone, sometimes even fractures may occur
- Damage to the lateral ligaments
- Incorrect position of the components
- Incorrect resection

#### After surgery

- Cardiovascular defects, such as TEN, deep phlebothrombosis, MI, postoperative hematoma
- Wound healing disorders, infection
- Loosening of the entire total replacement or its components, deformation or breakage of a component, flexion contracture, limb shortening

### Product disposal

After being used, the implant is classified as "dangerous waste" in health care. It is to be disposed of in accordance with the valid legislation of the country where the implant is being disposed of.

### Manufacturer's recommendations:

The implant cannot be re-sterilised. If the sterile packaging is damaged, do not use the medical device and send it back to the manufacturer.

### Risks arising from repeated use from the patient's perspective:

- Failure or instability of the implants due to damage to the articulated and anchoring surfaces (risk of rubbing, effect of rubbed parts and subsequent loosening of the implant) or inadequate sterilization, pain, risk of premature need to reoperate
- Risk of infection in the patient (risk of infection and loosening of the implant) – risk of urgent need to reoperate
- Mistaken implant size or type or incomplete implant – risk of repeated surgery or premature need to reoperate

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