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Mail: mailbox@beznoska.cz

<http://www.beznoska.cz>

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Procedure for cleaning, disinfection, inspection, sterilization and storage of non-sterile instruments and kits

1/ General information

- Dirty, unsterile or damaged instruments must not be used under any circumstances.
- Each instrument must be inspected before use, if any damage is found, it cannot be used.
- Instruments must not be modified or otherwise processed.
- When handling an instrument that has sharp edges or points, extra care must be taken. When using protective gloves there is a risk of damage.
- After cleaning/disinfecting/sterilizing, it is necessary to check that the cutting surfaces, rasps and other cutting tools are sufficiently sharp.
- Instruments should be checked for bone/tissue or body fluid residues after cleaning. If residues are evident, cleaning should be repeated.
- Instruments must not be exposed to substances containing chlorine or fluorine.
- Instruments must not be placed in saline solutions as corrosion and surface changes can be expected with prolonged exposure.
- Instruments made of titanium or titanium alloys are susceptible to discoloration caused by impurities from water vapors and cleaning agent residues, which form multi-colored oxide deposit layers on the surface. These oxide layers can darken after repeated sterilization and impair the legibility of the order number, LOT or batch (see Chapter 3. Materials). The washing and sterilization procedure described in these instructions must be followed.
- All surfaces must be accessible for sterilizing agents.
- Flash sterilisation is only possible for unpackaged medical devices intended for immediate use.
- All serious defects occurring in connection with instruments supplied by BEZNOSKA, s.r.o. must be reported immediately to mailbox@beznoska.cz.
- BEZNOSKA, s.r.o. instruments must be disposed of in accordance with national and international regulations.

- The cleaning, disinfection, inspection, sterilization and storage of instruments is the responsibility of the provider, i.e. the hospital.

2/ How not to proceed

Any residue that sticks and dries on the instruments leads to corrosion. Exposing instruments to moisture for extended periods of time also leads to damage!

Possible initial and subsequent defects and their causes:

Cause	Defect
Blood, pus, secretions, tissue debris, bone debris	Corrosion, rusting
Saline solution, iodine tincture, unsuitable water, unsuitable and/or improperly used cleaning agents and disinfectants	Pitting corrosion, discolouration
Steel wire and/or steel brush damage	Contact corrosion, material surface destruction, oxide layer removal, leading to increased susceptibility
Physical contact between instruments made of different metals	Contact corrosion
Instrument overload	Cutting surfaces became dull, damaged → increased vulnerability to corrosion
Physical contact between any instruments	Cutting surfaces become dull, damage and increased susceptibility to corrosion
Dirt in the sterilizer, e.g. due to corroded instruments or improper maintenance of the sterilizer.	Starting rust: contamination of intact tools with rust
Insufficient drying of instruments	Corrosion, rusting

How to avoid bigger issues:

- Never let residue (blood, secretions, tissue debris) dry on the instrument; clean it immediately after the surgery.
- Disassemble the instruments, put special care for cleaning of the recesses, use soft brushes.
- Use only cleaning agents and disinfectants intended for the specific materials and follow the manufacturer's instructions for use.
- Very carefully rinse disinfectants and cleaning agents with water.

- Never leave or store instruments when wet or damp.

3/ Materials

Stainless steel

Stainless steel corrosion resistance is due to the formation of a passive layer (a layer of chromium oxide) on its surface. The passive layer is extremely resistant to many chemical materials and physical parameters. However, it is a mistake to think that "stainless" steel cannot corrode. This material can be affected by external conditions such as lack of or improper care. The use of disinfectants and cleaning agents containing one or more of the following is not recommended for stainless steel: chlorine, oxalic acid, hydrogen peroxide (H₂O₂). If this is not followed, pitting and contact corrosion may occur, or the material may discolour.

Wrought titanium alloy

Wrought titanium alloy is a material that is very resistant to corrosion and external conditions due to the self-oxidation on its surface. The use of disinfectants and cleaning agents containing one or more of the following is not recommended for titanium alloy: chlorine, oxygenated acids (e.g. nitric acid, sulfuric acid, oxalic acid), hydrogen peroxide (H₂O₂). If this is not followed, the material may discolour.

Plastics

The plastics used for BEZNOSKA, s.r.o. instruments are very durable and can be sterilized at temperatures up to 134 °C (273 °F). The use of disinfectants and cleaning agents containing one or more of the following ingredients is not recommended for plastics: organic solvents (alcohols, ethers, ketones and benzenes), hydrogen peroxide (H₂O₂), aldehyde, halogens (chlorine, iodine, bromine). If this is not followed, the plastic may be deformed and destroyed.

Important notice:

Instruments must never be cleaned with a metal or steel brush or wire. None of the instruments must be exposed to temperatures above 134°C.

4/ Responsibility of the workplace – hospital

- Provides highly sterile surgical conditions.
- Ensures sterility of used instruments.
- Ensures the completeness and functionality of instruments for implantation.
- Ensures the usage of instruments supplied by the manufacturer .
- Ensures that the instruments are used for the purpose for which it is intended.
- Ensures that instruments are not modified or altered.
- Ensures that cleaning/disinfection/sterilization equipment used is regularly serviced and inspected:
 - Equipment should comply with ISO 15883 and bear the CE marking.
 - Use a program that is appropriate for the instruments.
 - Use sterile water, low microbial water (max. 10microbes/ml) or low endotoxin water (max. 0.25 endotoxin units/ml - purified water, highly purified water).
 - For drying use filtered air.
- Ensures that a suitable washing detergent for metal and plastic instruments is used (see chapter 3. Materials).
- Ensures that generally accepted safety precautions are followed when handling contaminated or potentially contaminated tools, such as the use of protective equipment - clothing, gloves, slippers, goggles, respirator.
- Ensures compliance with national hygiene legislation in designated country.
- Ensures monitoring and recording of sterilization parameter settings.

5/ Pre-sterilisation preparation

Pre-sterilization preparation is divided into manual preparation and automatic cleaning.

5.1/ Manual preparation

Rough contamination must be removed from the instruments immediately after use (within two (2) hours at the latest).

Multi-part instruments must first be disassembled into their individual parts. Instruments that are channelled (have a cavity) must first be cleaned with soft brushes (use the brushes that are designed for this purpose, note the warning in point 3). The brush must be inserted into the cavity at least 5 times, followed by flushing the cavity with water (see point 4 for water parameters) from a pressure gun, followed by blowing the cavity with compressed clean air. The moving parts of the instrument should be moved back and forth several times to wash the entire range of movement. The instruments can then be placed in a water bath with an enzyme solution or in a disinfectant solution for 10 min. In case of heavy contamination, the solutions should be repeatedly changed. It is possible to clean the instruments in an ultrasonic bath at a frequency of 35 - 40kHz for 5 min. It is necessary to proceed in such a way that the instruments are not damaged. At the end of the manual cleaning, the instruments should be thoroughly rinsed with clean water, repeating this process at least twice (see point 4 for water parameters).

5.2/ Automatic cleaning

The washing and disinfection machine is set to a standard cycle consisting of pre-wash, cleaning (washing), several rinses, followed by thermal disinfection and drying. Use a detergent that meets the material requirements (see point 3. Materials). Improperly selected detergent may cause staining, corrosion or discoloration of the base material of the instrument. The tool trays are washable and conform to ISO 11607.

5.3/ Warning after cleaning

Cleaned and disinfected instruments must be kept dry in the inter-operative period (moisture can cause corrosion of the instruments). Blunt, damaged and corroded instruments must be discarded; further use is not permitted. Corroded instruments must not be sterilised (risk of contamination of the autoclave and other instruments with rust particles and subsequent spread of corrosion). Before sterilization, check all instruments to ensure that all contaminants have been completely removed from the surface. If dirt is still visible, repeat the cleaning process again.

6/ Sterilization

Only the sterilization procedures listed below can be used to sterilize reusable surgical instruments. Other sterilization procedures are prohibited!

Steam sterilization:

- Fractionation vacuum method or gravity method (with sufficient drying of the equipment).
- ISO 15883 compliant steam sterilizer.
- The trays in which the instruments are supplied are also suitable for the sterilization cycle and comply with ISO 11607.

Methods of sterilization:

- Sterilization temperature 121 °C
- Overpressure of 105kPa (1,05 bar)
- Pressure of 205kPa (2,05bar)
- Sterilization exposure time 20 min

Or

Flash sterilization (for unpacked instruments only)

- Sterilization temperature 134 °C
- Overpressure of 204kPa (2,04 bar)
- Pressure of 304kPa (3,04 bar)
- Sterilization exposure time 4 min

Or

- Sterilization temperature 134 °C
- Overpressure of 204kPa (2,04 bar)
- Pressure of 304kPa (3,04 bar)
- Sterilization exposure time 7 min

Or

- Sterilization temperature 134 °C
- Overpressure of 204kPa (2,04 bar)
- Pressure of 304kPa (3,04 bar)
- Sterilization exposure time 3 min

Subsequently, provided that the temperature, overpressure, and pressure are maintained, the sterilization duration can be selected according to national regulations (legislation). However, the minimum exposure time is 3 minutes at the above-mentioned parameters to ensure the sterility of reusable surgical instruments.

Important notice:

Do not use hot air sterilization, radiation sterilization, formaldehyde sterilization, ethylene oxide sterilization or plasma sterilization.

Flash sterilization must not be used to sterilize instruments with a cavity.

Always follow your sterilizer manufacturer's instructions for use, especially with regard to cartridge weight, operating time and functional testing.

Corroded and rusted instruments can contaminate the water circulation in the sterilizer with rust particles. These rust particles will cause starting corrosion on intact instruments for all future sterilization cycles. It is therefore important to check and clean the unit regularly!

7/ Storage and transport of sterilised material

Packages with sterilised instrument/repair shall be stored either loose with a short expiration time or with a longer expiration time in a sterility-proof container that protects against dust in a closed cabinet, storage container, drawer or other container. For storage, a temperature of 15 to 30 °C and a relative humidity of 30 to 70 % is specified.

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